


RX Biliary System™

CANNULAS







Injection Lumens



Tandem™ RX Cannula

Order Number	Description	Shaft O.D. (F) / (mm)	Length (cm)	Recommended Guidewire (in) / (mm)
M00545820	Tandem RX Cannula	7-5.5 / 2.33-1.84	210	.035 / 0.89

Packaged one per box



RX Cannulas

Order Number	Description	Shaft O.D. (F) / (mm)	Length (cm)	Recommended Guidewire (in) / (mm)
M00545760	RX Standard ERCP Cannula	8.5-5 / 2.83-1.67	210	.035 / 0.89
M00545770	RX Tapered ERCP Cannula	8.5-5 / 2.83-1.67	210	.035 / 0.89
M00545780	RX Ball Tip ERCP Cannula	8.5-5 / 2.83-1.67	210	.035 / 0.89

Packaged one per box

RX Biliary System™

SPHINCTEROTOMES



Dreamtome™ RX Cannulating Sphincterotomes with Pre-loaded Dreamwire™ Guidewire




Order Number	Description	Cut Wire (mm)	Tip Length (mm)	Tip Diameter (F) / (mm)	Pre-loaded Guidewire** Length (cm)	O.D. (in) / (mm)
M00584000	Dreamtome RX 49	20	5	4.9 / 1.63	260	.035 / 0.89
M00584010	Dreamtome RX 49	30	5	4.9 / 1.63	260	.035 / 0.89
M00584020	Dreamtome RX 49	20	5	4.9 / 1.63	450	.035 / 0.89
M00584030	Dreamtome RX 49	30	5	4.9 / 1.63	450	.035 / 0.89
M00584040*	Dreamtome RX 44	20	5	4.4 / 1.47	260	.035 / 0.89
M00584050*	Dreamtome RX 44	30	5	4.4 / 1.47	260	.035 / 0.89
M00584060*	Dreamtome RX 44	20	5	4.4 / 1.47	450	.035 / 0.89
M00584070*	Dreamtome RX 44	30	5	4.4 / 1.47	450	.035 / 0.89

Packaged one per box

* Includes handle brake.

** Pre-loaded with Straight-tip guidewires only.



Hydratome™ RX Cannulating Sphincterotomes with Pre-loaded Hydra Jagwire™ Guidewire



Order Number	Description	Cut Wire (mm)	Shaft O.D. (F) / (mm)	Tip Length (mm)	Tip Diameter (F) / (mm)	Catheter Length (cm)	Pre-loaded Guidewire** Length (cm)	O.D. (in) / (mm)
M00583000	Hydratome RX 49	20	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	260	.035 / 0.89
M00583010	Hydratome RX 49	30	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	260	.035 / 0.89
M00583020	Hydratome RX 49	20	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	450	.035 / 0.89
M00583030	Hydratome RX 49	30	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	450	.035 / 0.89
M00583040*	Hydratome RX 44	20	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	260	.035 / 0.89
M00583050*	Hydratome RX 44	30	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	260	.035 / 0.89
M00583060*	Hydratome RX 44	20	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	450	.035 / 0.89
M00583070*	Hydratome RX 44	30	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	450	.035 / 0.89

Packaged one per box

* Includes handle brake.

** Pre-loaded with .035in (0.89mm) Hydra Jagwire™ Guidewire.

RX Biliary System™

SPHINCTEROTOMES (CONTINUED)



Jagtome™ RX Cannulating Sphincterotomes with Pre-loaded Jagwire™ Guidewire



Order Number	Description	Cut Wire (mm)	Shaft O.D. (F) / (mm)	Tip Length (mm)	Tip Diameter (F) / (mm)	Catheter Length (cm)	Pre-loaded Guidewire** Shaft Stiffness	Pre-loaded Guidewire** Length (cm)	O.D. (in) / (mm)
M00573000	Jagtome RX 49	20	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	Standard	260	.035 / 0.89
M00573020	Jagtome RX 49	30	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	Standard	260	.035 / 0.89
M00573010	Jagtome RX 49	20	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	Standard	450	.035 / 0.89
M00573030	Jagtome RX 49	30	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	Standard	450	.035 / 0.89
M00573500	Jagtome RX 49	20	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	Stiff	260	.035 / 0.89
M00573520	Jagtome RX 49	30	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	Stiff	260	.035 / 0.89
M00573510	Jagtome RX 49	20	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	Stiff	450	.035 / 0.89
M00573530	Jagtome RX 49	30	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	Stiff	450	.035 / 0.89
M00573040*	Jagtome RX 44	20	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	Standard	260	.035 / 0.89
M00573060*	Jagtome RX 44	30	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	Standard	260	.035 / 0.89
M00573050*	Jagtome RX 44	20	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	Standard	450	.035 / 0.89
M00573070*	Jagtome RX 44	30	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	Standard	450	.035 / 0.89
M00573540*	Jagtome RX 44	20	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	Stiff	260	.035 / 0.89
M00573560*	Jagtome RX 44	30	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	Stiff	260	.035 / 0.89
M00573550*	Jagtome RX 44	20	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	Stiff	450	.035 / 0.89
M00573570*	Jagtome RX 44	30	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	Stiff	450	.035 / 0.89
M00573580*	Jagtome RX 39	20	7-5.5 / 2.33-1.83	5	3.9 / 1.30	200	Standard	260	.025 / 0.64
M00573600*	Jagtome RX 39	30	7-5.5 / 2.33-1.83	5	3.9 / 1.30	200	Standard	260	.025 / 0.64
M00573590*	Jagtome RX 39	20	7-5.5 / 2.33-1.83	5	3.9 / 1.30	200	Standard	450	.025 / 0.64
M00573610*	Jagtome RX 39	30	7-5.5 / 2.33-1.83	5	3.9 / 1.30	200	Standard	450	.025 / 0.64
M00573080*	Jagtome RX 39	20	7-5.5 / 2.33-1.83	5	3.9 / 1.30	200	Stiff	260	.025 / 0.64
M00573100*	Jagtome RX 39	30	7-5.5 / 2.33-1.83	5	3.9 / 1.30	200	Stiff	260	.025 / 0.64



Packaged one per box

* Includes handle brake.

** Pre-loaded with the high performance Jagwire™ Guidewire (Jagtome RX 49 and RX 44 contain the Stiff Shaft Jagwire and Jagtome RX 39 the Standard Jagwire).

RX Biliary System™

SPHINCTEROTOMES (CONTINUED)



Autotome™ RX Cannulating Sphincterotomes

Order Number	Description	Cut Wire (mm)	Shaft O.D. (F) / (mm)	Tip Length (mm)	Tip Diameter (F) / (mm)	Catheter Length (cm)	Recommended Guidewire (in) / (mm)
M00545150	Autotome RX 49	20	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	.035 / 0.89
M00545160	Autotome RX 49	30	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	.035 / 0.89
M00545170*	Autotome RX 44	20	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	.035 / 0.89
M00545180*	Autotome RX 44	30	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	.035 / 0.89
M00545190*	Autotome RX 39	20	7-5.5 / 2.33-1.83	5	3.9 / 1.30	200	.025 / 0.64
M00545200*	Autotome RX 39	30	7-5.5 / 2.33-1.83	5	3.9 / 1.30	200	.025 / 0.64

Packaged one per box

Recommended Guidewires: For Autotome RX 44 and 49, recommended guidewire is .035in (0.89mm) Dreamwire™ High Performance Guidewire
For Autotome RX 39, recommended guidewire is .025in (0.69mm) Jagwire™ Guidewire

* Includes handle brake.

RX Triple-lumen Needleknife


Order Number	Description	Cut Wire (mm)	Shaft O.D. (F) / (mm)	Tip Length (mm)	Tip Diameter (F) / (mm)
M00545840	Needleknife RX	4-6	7-5.5 / 2.3-1.83	5	5 / 1.67

Packaged one per box

Recommended Guidewire: .035 in (0.89mm) Dreamwire™ High Performance Guidewire

RX Biliary System™

STONE RETRIEVAL




NEW Extractor™ Pro RX Triple-lumen Retrieval Balloons – Distal Injection



Order Number	Description	Balloon O.D. (mm)	Injection	Catheter O.D. (F) / (mm)	System Length (cm)
M00547000	Extractor Pro RX Triple-lumen Retrieval Balloon	9.0-12.0	Above	7-6 / 2.33-2.00	200
M00547010	Extractor Pro RX Triple-lumen Retrieval Balloon	12.0-15.0	Above	7-6 / 2.33-2.00	200
M00547020	Extractor Pro RX Triple-lumen Retrieval Balloon	15.0-18.0	Above	7-6 / 2.33-2.00	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire




NEW Extractor™ Pro RX Triple-lumen Retrieval Balloons – Proximal Injection



Order Number	Description	Balloon O.D. (mm)	Injection	Catheter O.D. (F) / (mm)	System Length (cm)
M00547030	Extractor Pro RX Triple-lumen Retrieval Balloon	9.0-12.0	Below	7-6 / 2.33-2.00	200
M00547040	Extractor Pro RX Triple-lumen Retrieval Balloon	12.0-15.0	Below	7-6 / 2.33-2.00	200
M00547050	Extractor Pro RX Triple-lumen Retrieval Balloon	15.0-18.0	Below	7-6 / 2.33-2.00	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire



Trapezoid™ RX Wire-guided Retrieval Baskets



Order Number	Description	Open Basket Diameter (cm)	Required Working Channel (mm)
M00510860	Trapezoid RX Wire-guided Retrieval Basket	1.5	3.2
M00510870	Trapezoid RX Wire-guided Retrieval Basket	2.0	3.2
M00510880	Trapezoid RX Wire-guided Retrieval Basket	2.5	3.2
M00510890*	Trapezoid RX Wire-guided Retrieval Basket	3.0	3.2
M00550620	Alliance™ II Inflation Handle**		

Packaged one per box


Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

* The Trapezoid 3 x 6cm basket is designed for crushing calculus larger than 1.5cm in diameter.

** Alliance II recommended for use as Mechanical Lithotripter. Alliance II is reusable.

RX Biliary System™

DILATION



Hurricane™ RX Biliary Dilation Balloons

Order Number	Description	Balloon O.D. (mm)	Balloon Length (cm)	Catheter O.D. (F) / (mm)	Maximum Inflation Pressure (ATM) / (kPa)	Minimum Channel Size (mm)	Recommended Guidewire (in) / (mm)
M00545890	Hurricane RX Dilation Balloon	4	2	5.8 / 1.93	12 / 1216	3.2	.035 / 0.89
M00545900	Hurricane RX Dilation Balloon	4	4	5.8 / 1.93	12 / 1216	3.2	.035 / 0.89
M00545910	Hurricane RX Dilation Balloon	6	2	5.8 / 1.93	11 / 1114	3.2	.035 / 0.89
M00545920	Hurricane RX Dilation Balloon	6	4	5.8 / 1.93	11 / 1114	3.2	.035 / 0.89
M00545930	Hurricane RX Dilation Balloon	8	2	5.8 / 1.93	11 / 1114	3.2	.035 / 0.89
M00545940	Hurricane RX Dilation Balloon	8	4	5.8 / 1.93	11 / 1114	3.2	.035 / 0.89
M00545950	Hurricane RX Dilation Balloon	10	2	5.8 / 1.93	8 / 810	3.2	.035 / 0.89
M00545960	Hurricane RX Dilation Balloon	10	4	5.8 / 1.93	8 / 810	3.2	.035 / 0.89

Packaged one per box

Recommended Inflation Device: Encore™ 26 Inflation Device, M00566670

CYTOLOGY



Combo Cath™ RX Cytology Brush

Order Number	Description	Brush Width (mm)	Catheter O.D. (F) / (mm)	Working Length (cm)	Recommended Guidewire (in) / (mm)
M00545000	Combo Cath RX Cytology Brush	2.1	8 / 2.67	200	.035 / 0.89

Packaged one per box


ACCESSORIES

RX Accessories



Order Number	Description	Units
M00545471	RX Locking Device	box 10
M00545211†	RX Locking Device – Olympus/Fujinon Brand Endoscope Compatible	box 10
M00545221†	RX Locking Device – Pentax Brand Endoscope Compatible	box 10
M00545261	RX Locking Device and Biopsy Cap System – Olympus Compatible	box 10
M00545271	RX Locking Device and Biopsy Cap System – Pentax Compatible	box 10
M00545281	RX Locking Device and Biopsy Cap System – Fujinon Compatible	box 10
M00566001	Rapid Refill Continuous Injection System	box 10

† Product not available in the European Economic Area (EEA), except in Greece. Not for use or distribution in the EEA, except in Greece. Please check availability of the product with your local sales representative or your local customer service.



Encore™ 26 Inflation Device




Order Number	Description	Units
M00566670	Encore 26 Inflation Device	box 1

RX Biliary Metal and Plastic Stents are located in the Stents Section on pages 19–20 and 22–25

Biliary Guidewires are located in the Guidewires Section on page 47–48

RX Biliary System – While in Stock

STONE RETRIEVAL




Extractor™ RX Retrieval Balloons – Distal Injection



Order Number	Description	Balloon O.D. (mm)	Injection	Catheter O.D. (F) / (mm)
M00546900	Extractor RX Retrieval Balloon	9.0-12.0	Above	7-6 / 2.33-2.00
M00546910	Extractor RX Retrieval Balloon	12.0-15.0	Above	7-6 / 2.33-2.00
M00546920	Extractor RX Retrieval Balloon	15.0-18.0	Above	7-6 / 2.33-2.00

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire



Extractor™ RX Retrieval Balloons – Proximal Injection




Order Number	Description	Balloon O.D. (mm)	Injection	Catheter O.D. (F) / (mm)
M00546930	Extractor RX Retrieval Balloon	9.0-12.0	Below	7-6 / 2.33-2.00
M00546940	Extractor RX Retrieval Balloon	12.0-15.0	Below	7-6 / 2.33-2.00
M00546950	Extractor RX Retrieval Balloon	15.0-18.0	Below	7-6 / 2.33-2.00

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

Biliary

CANNULAS




Tandem™ XL Triple-lumen ERCP Cannula



Order Number	Description	Shaft O.D. (F) / (mm)	Distal O.D. (F) / (mm)	Length (cm)	Tip	Recommended Guidewire (in) / (mm)
M00535700	Tandem XL Triple-lumen ERCP Cannula	7.0-5.5 / 2.33-1.83	5.5 / 1.83	210	Tapered	.035 / 0.89

Packaged one per box



Contour™ ERCP Cannulas




Order Number	Description	Shaft O.D. (F) / (mm)	Length (cm)	Tip	Recommended Guidewire (in) / (mm)
M00530850	Contour ERCP Cannula	5 / 1.67	210	Standard	.035 / 0.89
M00530860	Contour ERCP Cannula	5 / 1.67	210	Tapered	.035 / 0.89
M00530870	Contour ERCP Cannula	5 / 1.67	210	Ultra Tapered	.025 / 0.64
M00530880	Contour ERCP Cannula	5 / 1.67	210	5-4-3	.018 / 0.46
M00530890	Contour ERCP Cannula	5 / 1.67	210	Ball Tip	.035 / 0.89

Packaged one per box

Biliary

CANNULAS (CONTINUED)




Fluoro Tip™ ERCP Cannulas



Order Number	Description	Shaft O.D. (F) / (mm)	Length (cm)	Recommended Guidewire (in) / (mm)
M00530960	Fluoro Tip ERCP Cannula, Standard Tip	5.5 / 1.83	210	.035 / 0.89
M00530970	Fluoro Tip ERCP Cannula, Tapered Tip	5.5 / 1.83	210	.035 / 0.89
M00530980	Fluoro Tip ERCP Cannula, Ultra Tapered Tip	5.5 / 1.83	210	.025 / 0.64

Packaged one per box

SPHINCTEROTOMES



Stonetome™ Stone Removal Device




Order Number	Description	Cut Wire (mm)	Balloon Diameter (mm)	Tip Length (mm)	Catheter Length (cm)
M00535150	Stonetome (below cut wire)	20	11.5	5	200
M00535170	Stonetome (below cut wire)	20	11.5	20	200
M00535190	Stonetome (below cut wire)	30	11.5	5	200
M00535210	Stonetome (below cut wire)	30	11.5	20	200
M00535110	Stonetome (above cut wire)	20	11.5	20	200
M00535130	Stonetome (above cut wire)	30	11.5	20	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

Biliary

SPHINCTEROTOMES (CONTINUED)




Ultratome XL™ Triple-lumen Wire-guided Sphincterotomes

Order Number	Description	Cut Wire (mm)	Distal O.D. (F) / (mm)	Tip Length (mm)	Catheter Length (cm)
M00535900	Ultratome XL Short Nose	20	5.5 / 1.83	5	200
M00535910	Ultratome XL Long Nose	20	5.5 / 1.83	20	200
M00535920	Ultratome XL Short Nose	30	5.5 / 1.83	5	200
M00535930	Ultratome XL Long Nose	30	5.5 / 1.83	20	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire




Ultratome™ Double-lumen Wire-guided Sphincterotomes

Order Number	Description	Cut Wire (mm)	Distal O.D. (F) / (mm)	Tip Length (mm)	Catheter Length (cm)
M00530800	Ultratome Short Nose	20	5.5 / 1.83	5	200
M00530810	Ultratome Long Nose	20	5.5 / 1.83	20	200
M00530820	Ultratome Short Nose	30	5.5 / 1.83	5	200
M00530830	Ultratome Long Nose	30	5.5 / 1.83	20	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire



Microknife™ XL Triple-lumen Needleknife


Order Number	Description	Shaft O.D. (F) / (mm)	Catheter Length (cm)
M00532810	Microknife XL Triple-lumen Needleknife	7-5 / 2.33-1.83	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

Biliary

STONE RETRIEVAL




NEW Extractor™ Pro XL Triple-lumen Retrieval Balloons – Distal Injection

Order Number	Description	Balloon O.D. (mm)	Injection	Catheter O.D. (F) / (mm)	System Length (cm)
M00547100	Extractor Pro XL Triple-lumen Retrieval Balloon	9.0-12.0	Above	7-6 / 2.33-2.00	200
M00547110	Extractor Pro XL Triple-lumen Retrieval Balloon	12.0-15.0	Above	7-6 / 2.33-2.00	200
M00547120	Extractor Pro XL Triple-lumen Retrieval Balloon	15.0-18.0	Above	7-6 / 2.33-2.00	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire



NEW Extractor™ Pro XL Triple-lumen Retrieval Balloons – Proximal Injection

Order Number	Description	Balloon O.D. (mm)	Injection	Catheter O.D. (F) / (mm)	System Length (cm)
M00547130	Extractor Pro XL Triple-lumen Retrieval Balloon	9.0-12.0	Below	7-6 / 2.33-2.00	200
M00547140	Extractor Pro XL Triple-lumen Retrieval Balloon	12.0-15.0	Below	7-6 / 2.33-2.00	200
M00547150	Extractor Pro XL Triple-lumen Retrieval Balloon	15.0-18.0	Below	7-6 / 2.33-2.00	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

NEW Extractor™ Pro DL Double-lumen Retrieval Balloons


Order Number	Description	Balloon O.D. (mm)	Injection	Catheter O.D. (F) / (mm)	System Length (cm)
M00547200	Extractor™ Pro DL Double-lumen Retrieval Balloon	9.0-12.0	N/A	7-6 / 2.33-2.00	200
M00547210	Extractor™ Pro DL Double-lumen Retrieval Balloon	12.0-15.0	N/A	7-6 / 2.33-2.00	200
M00547220	Extractor™ Pro DL Double-lumen Retrieval Balloon	15.0-18.0	N/A	7-6 / 2.33-2.00	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

Biliary

DILATION



Max Force™ Biliary Dilation Balloons

Order Number	Description	Inflated O.D. (mm)	Balloon Length (cm)	Catheter Size (F) / (mm)	Required Working Channel (mm)
M00567330	Max Force Dilation Balloon	4	2	5 / 1.67	2.8
M00567340	Max Force Dilation Balloon	4	4	5 / 1.67	2.8
M00567370	Max Force Dilation Balloon	6	2	5 / 1.67	2.8
M00567380	Max Force Dilation Balloon	6	4	5 / 1.67	2.8
M00567410	Max Force Dilation Balloon	8	2	5 / 1.67	2.8
M00567420	Max Force Dilation Balloon	8	3	5 / 1.67	2.8

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

Recommended Inflation Device: Encore™ 26 Inflation Device, M00566670

ACCESSORIES


ACCESSORIES

Order Number	Description	Units
M00561250	Active Cord, Fits Olympus Brand Generator	box 1
M00561270	Active Cord, Fits Bovie®, Endostat™, Endostat II, Valleylab, Aspen Labs, ERBE Brand Generators	box 1
M00566670	⊗ Encore 26 Inflation Device	box 1
M00566001	⊗ Rapid Refill Continuous Injection System	box 10

Biliary Metal and Plastic Stents are located in the Stents Section on pages 21 and 26–27

Biliary Guidewires are located in the Guidewires Section on page 47

STONE RETRIEVAL




Extractor™ XL Triple-lumen Retrieval Balloons – Distal Injection

Order Number	Description	Inflated O.D. (mm)	Catheter Size (F) / (mm)	Usable Length (cm)
M00550450	Extractor XL Triple-lumen Retrieval Balloon	8.5	7 / 2.33 taper to 5 / 1.67	210
M00550460	Extractor XL Triple-lumen Retrieval Balloon	11.5	7 / 2.33 taper to 5 / 1.67	210
M00550470	Extractor XL Triple-lumen Retrieval Balloon	15.0	7 / 2.33 taper to 5 / 1.67	210

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire



Extractor™ XL Triple-lumen Retrieval Balloons – Proximal Injection

Order Number	Description	Inflated O.D. (mm)	Catheter Size (F) / (mm)	Usable Length (cm)
M00551450	Extractor XL Triple-lumen Retrieval Balloon	8.5	7 / 2.33 taper to 5 / 1.67	210
M00551460	Extractor XL Triple-lumen Retrieval Balloon	11.5	7 / 2.33 taper to 5 / 1.67	210
M00551470	Extractor XL Triple-lumen Retrieval Balloon	15.0	7 / 2.33 taper to 5 / 1.67	210

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire



Double-lumen Retrieval Balloons

Order Number	Description	Inflated O.D. (mm)	Catheter Size (F) / (mm)	Usable Length (cm)
M00551860	Double-lumen Retrieval Balloon	11.5	7 / 2.33	210
M00551850	Double-lumen Retrieval Balloon	15.0	7 / 2.33	210

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire