

Endoscopic Ultrasound (EUS)

FINE NEEDLE ASPIRATION (FNA)



NEW Expect™ Endoscopic Ultrasound Aspiration Needle*



Order Number	Description	Needle Size (gauge) / (mm)	Minimum Working Channel (mm)	Sheath Diameter (mm)	Units
M00550000	Expect Endoscopic Ultrasound Aspiration Needle	19 / 1.10	2.8	1.83	box 1
M00550001	Expect Endoscopic Ultrasound Aspiration Needle	19 / 1.10	2.8	1.83	box 5
M00550010	Expect Endoscopic Ultrasound Aspiration Needle	22 / 0.72	2.4	1.65	box 1
M00550011	Expect Endoscopic Ultrasound Aspiration Needle	22 / 0.72	2.4	1.65	box 5
M00550020	Expect Endoscopic Ultrasound Aspiration Needle	25 / 0.52	2.4	1.52	box 1
M00550021	Expect Endoscopic Ultrasound Aspiration Needle	25 / 0.52	2.4	1.52	box 5

*Includes a 20cc syringe and one-way stopcock.

Working length: 137.5cm to 141.5cm, adjustable.

Needle length: 0cm to 8cm, adjustable.