


Polypectomy


SNARES



Rotatable Snares*

Order Number	Description	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00561821	Small Oval	Medium Stiff	13	2.4	195	2.8	box 10
M00561831	Medium Oval	Medium Stiff	20	2.4	195	2.8	box 10

*With Star Catheter Technology.




Sensation™ Short Throw Snares

Order Number	Description	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00562691	Small Oval	Flexible	13	2.4	240	2.8	box 5
M00562692	Small Oval	Flexible	13	2.4	240	2.8	box 20
M00562693	Small Oval	Flexible	13	2.4	240	2.8	box 40
M00560321	Small Oval	Medium Stiff	13	2.4	240	2.8	box 5
M00562671	Medium Oval	Flexible	27	2.4	240	2.8	box 5
M00562672	Medium Oval	Flexible	27	2.4	240	2.8	box 20
M00562673	Medium Oval	Flexible	27	2.4	240	2.8	box 40
M00560311	Medium Oval	Medium Stiff	27	2.4	240	2.8	box 5
M00562651	Large Oval	Medium Stiff	30	2.4	240	2.8	box 5
M00562652	Large Oval	Medium Stiff	30	2.4	240	2.8	box 20
M00560331	Medium Crescent	Medium Stiff	27	2.4	240	2.8	box 5



Polypectomy

SNARES (CONTINUED)



Captiflex™ Snares

Order Number	Description	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00562471	Extra Small Oval	Flexible	11	2.4	240	2.8	box 10
M00562421	Small Oval	Flexible	13	2.4	240	2.8	box 10
M00562422	Small Oval	Flexible	13	2.4	240	2.8	box 40
M00562401	Medium Oval	Flexible	27	2.4	240	2.8	box 10
M00562402	Medium Oval	Flexible	27	2.4	240	2.8	box 40




Captivator™ Snares

Order Number	Description	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00562301	Small Oval	Stiff	13	2.4	240	2.8	box 10
M00562321	Medium Oval	Stiff	27	2.4	240	2.8	box 10
M00562391	Large Oval (thin wire)	Flexible	30	2.4	240	2.8	box 10
M00561311	Large Oval	Medium Stiff	30	2.4	240	2.8	box 10
M00562451	Small Hexagonal	Stiff	13	2.4	240	2.8	box 10
M00562341	Medium Hexagonal	Stiff	27	2.4	240	2.8	box 10
M00562371	Medium Crescent	Stiff	27	2.4	240	2.8	box 10



Polypectomy


SNARES (CONTINUED)



Captivator™ II Snares



Order Number	Description	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00561291	Extra Large Rounded	Stiff	33	2.4	240	2.8	box 10




Profile™ Pediatric Snares



Order Number	Description	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00562531	Extra Small Oval	Flexible	11	1.9	240	2.0	box 10
M00562551	Small Oval	Flexible	13	1.9	240	2.0	box 10
M00562571	Medium Oval	Flexible	27	1.9	240	2.0	box 10

Polypectomy

HOT BIOPSY FORCEPS




NEW Radial Jaw™ 4 Hot Biopsy Forceps



Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00515031	Radial Jaw 4 Hot Biopsy Forceps	2.4	240	2.8	box 5
M00515032	Radial Jaw 4 Hot Biopsy Forceps	2.4	240	2.8	box 20
M00515033	Radial Jaw 4 Hot Biopsy Forceps	2.4	240	2.8	box 40

POLYP RETRIEVAL DEVICE



TWISTER™ Rotatable Polyp Retrieval Device




Order Number	Description	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00561401	TWISTER Rotatable Polyp Retrieval Device	13	2.5	230	2.8	box 5
M00561402	TWISTER Rotatable Polyp Retrieval Device	13	2.5	230	2.8	box 10
M00561411	TWISTER Rotatable Polyp Retrieval Device	22	2.5	230	2.8	box 5
M00561412	TWISTER Rotatable Polyp Retrieval Device	22	2.5	230	2.8	box 10
M00561421	TWISTER Rotatable Polyp Retrieval Device	25	2.5	230	2.8	box 5
M00561422	TWISTER Rotatable Polyp Retrieval Device	25	2.5	230	2.8	box 10
M00561431	TWISTER Rotatable Polyp Retrieval Device	32	2.5	230	2.8	box 5
M00561432	TWISTER Rotatable Polyp Retrieval Device	32	2.5	230	2.8	box 10

Active Cords

Order Number	Description	Units
M00561250	Active Cord, Fits Olympus® Brand Generator	box 1
M00561270	Active Cord, Fits Endostat™, Endostat™ II, Bovie®, ValleyLab®, Aspen Labs and ERBE® Brand Generators	box 1

HOT BIOPSY FORCEPS



Radial Jaw™ 3 Hot Biopsy Forceps

Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00515501	Radial Jaw 3 Hot Biopsy Forceps	2.2	240	2.8	box 5
M00515502	Radial Jaw 3 Hot Biopsy Forceps	2.2	240	2.8	box 20
M00515503	Radial Jaw 3 Hot Biopsy Forceps	2.2	240	2.8	box 40