



# Bronchial Thermoplasty (BT) System



## **NEW** Alair™ Catheter ATS 2-5

Order Number	Model	Description	Active Electrode Length (mm)	Tip Diameter (mm)	Minimum Working Channel (mm)
M005ATS25020	ATS 2-5	Alair BT Catheter – Global	5.0	~1.5	2.0
M005ATS25010	ATS 2-5	Alair BT Catheter – North America	5.0	~1.5	2.0

Packaged one per box



## **NEW** Alair™ RF Controller ATS 200

Order Number	Model	Description
M005ATS20000	ATS 200	Alair RF Controller – monopolar (return electrode not included)

Packaged one per box

## **NEW** Alair™ BT Controller Accessory Kits ATS 201\*


Order Number	Model	Description
M005ATS20110	ATS 201	North American Plug
M005ATS20120	ATS 201	British Type G
M005ATS20130	ATS 201	EU Hybrid Type E/F
M005ATS20140	ATS 201	Australian Type I
M005ATS20150	ATS 201	Danish Type K
M005ATS20160	ATS 201	Brazil IEC 60906-1
M005ATS20170	ATS 201	Chinese CPCS-CCC (Chinese 10 A/250 V)
M005ATS20180	ATS 201	Swiss Type J
M005ATS20190	ATS 201	Israel Type H

Packaged one per box

\* Including Power Cord with Molded Plug, Foot Switch, Operating Manual.

# Pulmonary


## BIOPSY FORCEPS



### **NEW** Radial Jaw™ 4 Pulmonary Biopsy Forceps




Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Units
M00515181	Radial Jaw 4 Standard Capacity	1.8	100	2.0	Blue	box 5
M00515182	Radial Jaw 4 Standard Capacity	1.8	100	2.0	Blue	box 20
M00515191	Radial Jaw 4 Standard Capacity w/Needle	1.8	100	2.0	Blue	box 5
M00515192	Radial Jaw 4 Standard Capacity w/Needle	1.8	100	2.0	Blue	box 20
M00515201	Radial Jaw 4 Large Capacity	2.4	100	2.8	Blue	box 5
M00515202	Radial Jaw 4 Large Capacity	2.4	100	2.8	Blue	box 20



### Radial Jaw™ 3 Pulmonary Biopsy Forceps




Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Units
M00515231	Radial Jaw 3	1.8	100	2.0	Blue	box 5
M00515241	Radial Jaw 3 w/Needle	1.8	100	2.0	Blue	box 5
M00515211	Radial Jaw 3 Large Capacity	2.2	100	2.8	Blue	box 5



### Radial Edge™ Pulmonary Biopsy Forceps



Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Units
M00515561	Radial Edge	1.8	100	2.0	Blue	box 5
M00515591	Radial Edge w/Needle	1.8	100	2.0	Blue	box 5




### Radial Jaw™ Pulmonary Biopsy Forceps



Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Units
M00512661	Radial Jaw	1.8	100	2.0	Blue	box 5
M00512691	Radial Jaw w/Needle	1.8	100	2.0	Blue	box 5

# Pulmonary

## TRANSBRONCHIAL ASPIRATION NEEDLES



### eXcelon™ Transbronchial Aspiration Needle with Syringe

Order Number	Description	Gauge	Needle		Needle Length (cm)	Catheter Length (cm)	Shaft O.D. (mm)	Units
			O.D. (mm)	I.D. (mm)				
M00564101	eXcelon Transbronchial Aspiration Needle w/Syringe	19	1.07	0.69	15	130	1.8	box 5
M00564111	eXcelon Transbronchial Aspiration Needle w/Syringe	20	0.90	0.58	15	130	1.8	box 5
M00564121	eXcelon Transbronchial Aspiration Needle w/Syringe	21	0.81	0.50	15	130	1.8	box 5

## CYTOLOGY BRUSHES




### Celebrity™ Endoscopic Cytology Brushes – Bronchoscope

Order Number	Description	Required Working Channel (mm)	Bristle O.D. (mm)	Sheath Length (cm)	Units
M00516001	Celebrity Cytology Brush	2.0	1.0	140	box 10
M00516011	Celebrity Cytology Brush	2.0	1.5	140	box 10
M00516071	Celebrity Cytology Brush	2.0	1.9	100	box 10
M00516151	Celebrity Cytology Brush	2.0	1.9	150	box 10

# Pulmonary

## AIRWAY RETRIEVAL BASKETS



### ZeroTip™ Airway Retrieval Baskets



Order Number	Description	Sheath O.D. (mm)	Working Length (cm)	Basket Opening (mm)	Sheath Materials
M00513200	ZeroTip Airway Retrieval Basket	0.8	120	12	Polyimide/PTFE
M00513210	ZeroTip Airway Retrieval Basket	1.0	120	16	Polyimide/PTFE

Packaged one per box

## PULMONARY DILATION BALLOONS



### CRE™ Pulmonary Dilation Balloons




Order Number	Description	Balloon Length (cm)	Initial Balloon O.D. (mm)	Intermediate Balloon O.D. (mm)	Maximum Balloon O.D. (mm)	Working Catheter Length (cm)
M00550300	CRE Balloon	5.5	12	13.5	15	75
M00550310	CRE Balloon	5.5	15	16.5	18	75
M00550320	CRE Balloon	5.5	18	19.0	20	75
M00550330	CRE Balloon	3.0	8	9.0	10	75
M00550340	CRE Balloon	3.0	10	11.0	12	75
M00550350	CRE Balloon	3.0	12	13.5	15	75

Packaged one per box

# Pulmonary

## INFLATION SYSTEM



### Alliance™ II Inflation System

Order Number	Description	Units
M00550620	Alliance™ II Inflation Handle	box 1
M00550601	⊗ 60ml Syringe/Gauge Assembly	box 5

**Airway Stents are located in the Stents Section on pages  
33–35 and 37–39**

**Pulmonary Guidewires are located in the Guidewires Section on page 48**