


# DILATATION PRODUCTS

## URETERAL BALLOON DILATATION CATHETERS



### ■ UroMax Ultra™ Balloon Dilatation Catheters

Order Number	Balloon Inflated O.D.		Balloon Length (cm)	Catheter Size		Usable Length (cm)	Recommended Guidewire	
	(Fr)	(mm)		(Fr)	(mm)		(in)	(mm)
M0062251000	12	4.0	4	5.8	1.9	75	.038	0.97
M0062251010	15	5.0	4	5.8	1.9	75	.038	0.97
M0062251020	18	6.0	4	5.8	1.9	75	.038	0.97
M0062251030	21	7.0	4	5.8	1.9	75	.038	0.97
M0062251040	24	8.0	4	5.8	1.9	75	.038	0.97
M0062251050	30	10.0	4	5.8	1.9	75	.038	0.97
M0062251060	12	4.0	6	5.8	1.9	75	.038	0.97
M0062251070	15	5.0	6	5.8	1.9	75	.038	0.97
M0062251080	18	6.0	6	5.8	1.9	75	.038	0.97
M0062251090	21	7.0	6	5.8	1.9	75	.038	0.97
M0062251100	24	8.0	8	5.8	1.9	75	.038	0.97
M0062251110	30	10.0	8	5.8	1.9	75	.038	0.97
M0062251150	12	4.0	10	5.8	1.9	75	.038	0.97
M0062251160	15	5.0	10	5.8	1.9	75	.038	0.97
M0062251170	18	6.0	10	5.8	1.9	75	.038	0.97
M0062251180	21	7.0	10	5.8	1.9	75	.038	0.97


### ■ UroMax Ultra™ Balloon Dilatation Catheters (Kits)

Order Number	Balloon Inflated O.D.		Balloon Length (cm)	Catheter Size		Usable Length (cm)	Recommended Guidewire	
	(Fr)	(mm)		(Fr)	(mm)		(in)	(mm)
M0062251200	12	4.0	4	5.8	1.9	75	.038	0.97
M0062251210	15	5.0	4	5.8	1.9	75	.038	0.97
M0062251220	18	6.0	4	5.8	1.9	75	.038	0.97
M0062251230	21	7.0	4	5.8	1.9	75	.038	0.97
M0062251240	24	8.0	4	5.8	1.9	75	.038	0.97
M0062251250	30	10.0	4	5.8	1.9	75	.038	0.97
M0062251260	12	4.0	6	5.8	1.9	75	.038	0.97
M0062251270	15	5.0	6	5.8	1.9	75	.038	0.97
M0062251280	18	6.0	6	5.8	1.9	75	.038	0.97
M0062251290	21	7.0	6	5.8	1.9	75	.038	0.97
M0062251300	24	8.0	8	5.8	1.9	75	.038	0.97
M0062251310	30	10.0	8	5.8	1.9	75	.038	0.97
M0062251350	12	4.0	10	5.8	1.9	75	.038	0.97
M0062251360	15	5.0	10	5.8	1.9	75	.038	0.97
M0062251370	18	6.0	10	5.8	1.9	75	.038	0.97
M0062251380	21	7.0	10	5.8	1.9	75	.038	0.97

Kit Includes: UroMax Ultra™ High Pressure Balloon Dilatation Catheter, Encore™ 26 Inflator

## DILATATION PRODUCTS

### URETERAL BALLOON DILATATION CATHETERS (Continued)



#### Passport™ Balloon Dilatation Catheter (Kits) contains LATEX


Order Number	Balloon Inflated O.D.		Balloon Length (cm)	Tip Length (cm)	Shaft Diameter		Usable Length (cm)	Minimum Working Channel	
	(Fr)	(mm)			(Fr)	(mm)		(Fr)	(mm)
M0062181200	12	4.0	2	2	3	1.0	150	3.4	1.1
M0062181210	12	4.0	2	4	3	1.0	150	3.4	1.1
M0062181300	12	4.0	4	2	3	1.0	150	3.4	1.1
M0062181310	12	4.0	4	4	3	1.0	150	3.4	1.1

Kit Includes: Passport Balloon Dilatation Catheter, Encore™ 26 Inflator

#### Passport™ Balloon Dilatation Catheter contains LATEX

Order Number	Balloon Inflated O.D.		Balloon Length (cm)	Tip Length (cm)	Shaft Diameter		Usable Length (cm)	Minimum Working Channel	
	(Fr)	(mm)			(Fr)	(mm)		(Fr)	(mm)
M0062181000	12	4.0	2	2	3	1.0	150	3.4	1.1
M0062181010	12	4.0	2	4	3	1.0	150	3.4	1.1
M0062181100	12	4.0	4	2	3	1.0	150	3.4	1.1
M0062181110	12	4.0	4	4	3	1.0	150	3.4	1.1

## INFLATION DEVICES



#### Encore™ 26 Inflator

Order Number	Description
M0067101131*	20ml Encore 26 Inflator with Pressure Gauge
M0067101140	20ml Encore 26 Inflator with Pressure Gauge

\*Packaged Five (5) Per Box

# DILATATION PRODUCTS

## URETERAL DILATORS


### ■ Nottingham™ One-Step Tapered Dilators with HydroPlus™ Coating

Order Number	Taper Size From/To (Fr) (mm)		Taper Length (cm)	Usable Length (cm)
M0062301000	6/10	2.0/3.3	4	70
M0062301010	6/12	2.0/4.0	4	70
M0062301100	10/16	3.3/5.3	4	70
M0062301150	12/18	4.0/6.0	4	70

### ■ Ureteral Dilator with HydroPlus™ Coating and Procedural Sheath Set

Order Number	Description
M0062501000	Ureteral Dilators and Procedural Sheath Set

Set Includes: One 6Fr/2.0mm, 8Fr/2.7mm, 10Fr/3.3mm, 12Fr/4.0mm and 14Fr/4.7mm C-Flex™ Radiopaque Dilator with HydroPlus™ Coating; One 16Fr/5.3mm O.D., 14Fr/4.7mm I.D. x 37.5cm Length Radiopaque PTFE Procedural Sheath; .038in/0.97mm PTFE Coated Guidewire



### ■ Ureteral Dilatation Set with HydroPlus™ Coating

Order Number	Description
M0062401000	Ureteral Dilatation Set

Set Includes: One 6Fr/2.0mm, 8Fr/2.7mm, 10Fr/3.3mm, 12Fr/4.0mm, 14Fr/4.7mm and 16Fr/5.3mm C-Flex Radiopaque Dilator with HydroPlus Coating and .038in/0.97mm PTFE Coated Guidewire

### ■ Ureteral Dilators with HydroPlus™ Coating

Order Number	O.D. Size (Fr) (mm)		Usable Length (cm)
M0062401100	6	2.0	70
M0062401110	8	2.7	70
M0062401120	10	3.3	70
M0062401130	12	4.0	70
M0062401140	14	4.7	70
M0062401150	16	5.3	70