


## GUIDEWIRES

### SENSOR™ GUIDEWIRES



#### ■ Sensor™ Dual-Flex PTFE-Nitinol Guidewires with Hydrophilic Tip (10cm Flexible Proximal End)

Order Number	Diameter		Length (cm)	Tip Design	
	(in)	(mm)		(Shape)	(Style)
M0066703081	.035	0.89	150	Straight	3cm Flexible
M0066703011	.035	0.89	150	Angled	3cm Flexible
M0066703121	.038	0.97	150	Straight	3cm Flexible
M0066703021	.038	0.97	150	Angled	3cm Flexible


Packaged Five (5) Per Box

#### ■ Sensor™ PTFE-Nitinol Guidewires with Hydrophilic Tip

Order Number	Diameter		Length (cm)	Tip Design	
	(in)	(mm)		(Shape)	(Style)
M0066703051	.035	0.89	150	Straight	3cm Flexible
M0066703061	.035	0.89	150	Angled	3cm Flexible
M0066703091	.038	0.97	150	Straight	3cm Flexible
M0066703101	.038	0.97	150	Angled	3cm Flexible

Packaged Five (5) Per Box

## HYDROPHILIC GUIDEWIRES




#### ■ ZIPwire™ Hydrophilic Guidewires

Order Number	Diameter		Length (cm)	Tip Shape	Shaft Stiffness
	(in)	(mm)			
M006630201B1	.025	0.64	150	Straight	Regular
M006630216B1	.025	0.64	150	Straight	Stiff
M006630205B1	.035	0.89	150	Straight	Regular
M006630206B1	.035	0.89	150	Angled	Regular
M006630222B1	.035	0.89	150	Straight	Stiff
M006630223B1	.035	0.89	150	Angled	Stiff
M006630208B1	.038	0.97	150	Straight	Regular
M006630209B1	.038	0.97	150	Angled	Regular
M006630225B1	.038	0.97	150	Straight	Stiff

Packaged Five (5) Per Box

# GUIDEWIRES

## ZEBRA™ UROLOGICAL GUIDEWIRES




### Zebra™ Urological Guidewires

Order Number	Diameter		Length (cm)	Tip Design	
	(in)	(mm)		(Shape)	(Style)
M0066701081	.025	0.64	150	Straight	3cm Flexible
M0066701091	.025	0.64	150	Angled	3cm Flexible
M0066701101	.032	0.81	150	Straight	3cm Flexible
M0066701121	.035	0.89	150	Straight	3cm Flexible
M0066701131	.035	0.89	150	Angled	3cm Flexible
M0066701141	.038	0.97	150	Straight	3cm Flexible
M0066701151	.038	0.97	150	Angled	3cm Flexible

Supplied with Torque Vise

Packaged Five (5) Per Box

## PTFE COATED GUIDEWIRES



### PTFE Coated Guidewires (Fixed Core)


Order Number	Diameter		Length (cm)	Tip Design	
	(in)	(mm)		(Shape)	(Style)
M0066201070	.025	0.64	150	Straight	3cm Flexible
M0066201200	.035	0.89	150	Straight	3cm Flexible
M0066201250	.035	0.89	150	Straight	15cm Flexible; Bentson-Type
M0066201000	.038	0.97	150	Straight	3cm Flexible
M0066201050	.038	0.97	150	Straight	15cm Flexible; Bentson-Type

### PTFE Coated Guidewires (Fixed Core, Stiff Bodied)

Order Number	Diameter		Length (cm)	Tip Design	
	(in)	(mm)		(Shape)	(Style)
M0066201080	.025	0.64	150	Straight	3cm Flexible
M0066201100	.035	0.89	150	Straight	3cm Flexible
M0066201110	.035	0.89	150	3mm "J"	
M0066201010	.038	0.97	150	Straight	3cm Flexible
M0066201020	.038	0.97	150	3mm "J"	

## GUIDEWIRES

### PTFE COATED GUIDEWIRES (Continued)



#### ■ PTFE Coated Guidewires (Moveable Core)

Order Number	Diameter		Length (cm)	Tip Design	
	(in)	(mm)		(Shape)	(Style)
M0066201300	.035	0.89	150	Straight	
M0066201210	.035	0.89	150	3mm "J"	
M0066201030	.038	0.97	150	Straight	
M0066201220	.038	0.97	150	3mm "J"	

#### ■ Amplatz Super Stiff™ PTFE Coated Guidewires

Order Number	Diameter		Length (cm)	Tip Design	
	(in)	(mm)		(Shape)	(Style)
M0066401041	.035	0.89	145	Straight	7cm Flexible, Bentson-Type
M0066401051	.035	0.89	145	3mm "J"	v182.00/bx
M0066401081	.035	0.89	145	Straight	3.5cm Flexible
M0066401061	.038	0.97	145	Straight	6cm Flexible, Bentson-Type
M0066401071	.038	0.97	145	3mm "J"	
M0066401091	.038	0.97	145	Straight	3.5cm Flexible

Packaged Five (5) Per Box

## LUBRIGLIDE™ COATED GUIDEWIRES

#### ■ Standard Lubriglide™ Coated Guidewires (Fixed Core)

Order Number	Diameter		Length (cm)	Tip Design	
	(in)	(mm)		(Shape)	(Style)
M0066001060	.032	0.81	150	Straight	18cm Coil with 3cm Flexible; Lunderquist-Type
M0066001200	.035	0.89	150	Straight	3cm Flexible
M0066001250	.035	0.89	150	Straight	3cm Flexible
M0066001050	.038	0.97	150	Straight	15cm Flexible; Bentson-Type

#### ■ Standard Lubriglide™ Coated Guidewires (Fixed Core, Stiff Bodied)


Order Number	Diameter		Length (cm)	Tip Design	
	(in)	(mm)		(Shape)	(Style)
M0066001020	.038	0.97	150	Straight	3cm Flexible

#### ■ Standard Lubriglide™ Coated Guidewires (Moveable Core)

Order Number	Diameter		Length (cm)	Tip Design	
	(in)	(mm)		(Shape)	(Style)
M0066001300	.035	0.89	150	Straight	
M0066001030	.038	0.97	150	Straight	

# GUIDEWIRES

## GUIDEWIRE ACCESSORIES



### ■ Guidewire Torque Vise

---

Order Number	Description
M0067301301	Torque Vise Steerable Handle

---

Packaged Twelve (12) Per Box