



URETERAL CATHETERS

TORQUEABLE CATHETERS



■ Dual Lumen™ Ureteral Catheter

Order Number	O.D. Size		Working Length (cm)	Recommended Guidewire		Injection Lumen	
	(Fr)	(mm)		(in)	(mm)	(in)	(mm)
M0064051000	10	3.3	50	.038	0.97	.050	1.7




■ Imager™ II Catheters

Order Number	Description	O.D. Size		Working Length (cm)
		(Fr)	(mm)	
M0064003001	C1	5	1.7	65
M0064003011	C2	5	1.7	65
M0064003021	Straight	5	1.7	65
M0064003031	Bern	5	1.7	65
M0064003041	JB1	5	1.7	65
M0064004021	Straight	5	1.7	100
M0064004031	Bern	5	1.7	100
M0064004041	JB1	5	1.7	100
M0064005031	Bern	5	1.7	40

Packaged Five (5) Per Box

URETERAL CATHETERS

URETERAL CATHETERS



■ Flexima™ Ureteral Catheters

Order Number	Description	Catheter O.D. Size		Working Length (cm)	Tip O.D. Size	
		(Fr)	(mm)		(Fr)	(mm)
M0064002001	Open End	4	1.3	70		
M0064002011	Open End	5	1.7	70		
M0064002021	Open End	6	2.0	70		
M0064002031	Open End	7	2.3	70		
M0064002041	Open End	8	2.7	70		
M0064002111	Cone Tip	5	1.7	70	8	2.7
M0064002121	Cone Tip	6	2.0	70	10	3.3
M0064002131	Cone Tip	7	2.3	70	12	4.0
M0064002141	Cone Tip	8	2.7	70	14	4.7
M0064002201	Whistle Tip	4	1.3	70		
M0064002211	Whistle Tip	5	1.7	70		
M0064002221	Whistle Tip	6	2.0	70		
M0064002301	Olive Tip	4	1.3	70		
M0064002311	Olive Tip	5	1.7	70		
M0064002361	Spiral Tip	5	1.7	70		
M0064002371	Spiral Tip	6	2.0	70		
M0064002411	Wedge Tip	5	1.7	70	14	4.7

Each Catheter Supplied with Injection Hub

Packaged Twenty (20) Per Box


■ Flexima™ Ureteral Catheters (Sets with .035in/0.89mm PTFE Bentson-Tip Guidewire)

Order Number	Description	Catheter O.D. Size		Working Length (cm)	Tip O.D. Size (Fr)
		(Fr)	(mm)		
M0064002510	Angled Tip Set	5	1.3	70	-
M0064002520	Angled Tip Set	6	1.7	70	-

Set Includes: Angled Tip Catheter, .035in/0.89mm PTFE Bentson-Tip Guidewire and Injection Hub

URETERAL CATHETERS


URETERAL CATHETERS (Continued)



Axxcess™ 6 French Ureteral Catheter

Order Number	Catheter O.D. Size (Fr) (mm)	Working Length (cm)
M0064001161	6 2.0	70


Packaged Ten (10) Per Box



C-Flex™ Flexible Tip Open End Ureteral Catheters (Sets)

Order Number	O.D. Size (Fr) (mm)	Working Length (cm)
M0064001400	6 2.0	70
M0064001410	7 2.3	70

Set Includes: C-Flex™ Flexible Tip Catheter, .038in/0.97mm Stainless Steel Guidewire and (2) Luer Lock Catheter Adaptors




C-Flex™ Firm Open End Ureteral Catheters (Sets)

Order Number	O.D. Size (Fr) (mm)	Working Length (cm)
M0064001500	4.8 1.2	70
M0064001510	6.0 2.0	70
M0064001520	7.0 2.3	70
M0064001530	8.0 2.7	70

Set Includes: Firm C-Flex™ Catheter, .038in/0.97mm Stainless Steel Guidewire and Two (2) Luer Lock Catheter Adaptors

URETERAL CATHETERS


URETERAL CATHETERS (Continued)



■ C-Flex™ Firm Stamey Open Tip Ureteral Catheters (Sets)

Order Number	O.D. Size (Fr) (mm)	Working Length (cm)
M0064001610	4.0 1.3	70
M0064001620	4.8 1.6	70
M0064001630	6.0 2.0	70
M0064001640	7.0 2.3	70

Set Includes: Firm C-Flex™ Catheter, Placement Stylet, and Two (2) Luer Lock Catheter Adaptors



■ C-Flex™ Ureteral Pigtail Catheters

Order Number	O.D. Size (Fr) (mm)	Increments of Side Ports (cm)	Working Length (cm)
M0064001700*	4.8 1.6	5	70
M0064001710	6.0 2.0	5	70
M0064001720	7.0 2.3	5	70
M0064001730	8.0 2.7	5	70

Set Includes: Firm C-Flex™ Pigtail Catheter, .038in/0.97mm Lubriglide™ Coated Guidewire and Two (2) Luer Lock Catheter Adaptors

*Set Includes: Firm C-Flex™ Pigtail Catheter, .038in/0.97mm PTFE Coated Guidewire and Two (2) Luer Lock Catheter Adaptors