


PROCEDURE ACCESSORIES

INFLATION DEVICES




■ Encore™ 26 Inflator

Order Number	Description
M0067101131*	20ml Encore 26 Inflator with Pressure Gauge
M0067101140	20ml Encore 26 Inflator with Pressure Gauge

*Packaged Five (5) Per Box

IRRIGATION, EVACUATION AND RESECTION ACCESSORIES



■ Urovac™ Bladder Evacuator


contains
LATEX

Order Number	Description
M0067301251	Urovac Bladder Evacuator with Adaptor

Packaged Ten (10) Per Box


PROCEDURE ACCESSORIES

IRRIGATION, EVACUATION AND RESECTION ACCESSORIES (Continued)



■ SAPS™ Single Action Pumping System

Order Number	Description
M0067201001	Single Action Pumping System
System Includes: 10ml Vacuum Syringe, Check Valve with Adjustable Length Connecting Tubes	
Packaged Five (5) Per Box	




■ SAPS™ CF Continuous Flow Single Action Pumping System

Order Number	Description
M0067201011	Single Action Pumping System CF
System Includes: 10ml Vacuum Syringe, Continuous Flow Adjustable Length Tubing System with Dual-Check Valve and Roller Clamp	

PROCEDURE ACCESSORIES

OTHER



■ UroLok™ II Scope Adaptor

Order Number	Description
M0067301401	UroLok II with Ball-fitting Adaptor
Packaged Ten (10) Per Box	

■ Gateway™ Advantage Y Adaptor

Order Number	Description
M0067301501	Gateway Advantage Y-Adaptor
Packaged Ten (10) Per Box	

■ Connecting Tubes for Drainage Bag

Order Number	Description
M0067401000	30cm Connecting Tube for 3 Fr/1 mm through 6Fr / 2.0mm Catheters
M0067401010	30cm Connecting Tube for 7 Fr/2.3 mm through 10Fr / 3.3mm Catheters
M0067401020	30cm "Y" Connecting Tube for Two 3 Fr/1 mm through 6Fr / 2.0mm Catheters

■ Ureteral Catheter Connectors

Order Number	Description
M0067401100	Connection for 3 Fr/1 mm through 6Fr / 2.0mm Catheters with Luer Lock Fitting
M0067401110	Connection for 7 Fr/2.3 mm through 10Fr / 3.3mm Catheters with Luer Lock Fitting

■ Drainage Bag Connectors

Order Number	Descriptions
M0067401200	Connection for Female Luer Lock Fitting to Drainage Bag
M0067401220	Connection for 3 Fr/1 mm through 6Fr / 2.0mm Ureteral Catheters to Drainage Bag